

Immunofluorescence staining (cells and tissues)

Introduction

Get started by giving your protocol a name and editing this introduction.

Materials

- › Antigen Retrieval Buffer
 - › 10mM Tri-Sodium Citrate (dihydrate)
 - › 5mM EDTA
 - › 0,05% Tween 20
 - › pH 6,0
- › Xylene
- › Ethanol
 - › 100% / 95% / 70% / 50%
- › PBS
 - › 140 mM NaCl
 - › 0,27 mM KCl
 - › 0,15 mM K_2HPO_4
 - › 7,1 mM $Na_2HPO_4 \times 2H_2O$
 - › pH 7,2-7,4
- › IF buffer
- › 10mM ethanolamine. in PBS
- › NABH₄ (1%)
- › Aurion BSA (1%) = Blocking Buffer
- › Donkey Serum/ Sheep Serum
- › Liquid Blocker Pen
- › Fluorescence mounting media

Procedure

1st Day: removing Parafin, rehydrate, ... and incubate 1st antibody

- ✓ 1. Rinse cells or tissue in PBS briefly to remove media components. **(cells only)**

- ✓ 2. Fix cells or tissue in 4% paraformaldehyde in PBS (freshly prepared) for 15 min at 4°C. **(cells only)**
- ✓ 3. Prepare 10% para. stock like this: 1g paraformaldehyde in 10 ml dH₂O Heat to ~60°C Add 1 drop of 1M NaOH and continue heating until most is dissolved. Filter stock thru 0.2 µm filter **(cells only)**
- ✓ 4. For removing the Parafin use Xylene for Rehydration use different concentrations of Ethanol procede the Following Protocol **(Tissues only)**
 - 15 min Xylene
 - 15 min Xylene
 - 15 min Xylene
 - 5 min Ethanol 100%
 - 5 min Ethanol 100%
 - 5 min Ethanol 95%
 - 5 min Ethanol 70%
 - 5 min Ethanol 50%
 - 5 min PBS
- ✓ 5. While Samples are in the PBS prepare Antigen retrieval Buffer and Heat the Presure cooker up to 120°C **(Tissues only)**

Boil the Samples under Pressure for 30 minutes and afterwards let them cool to room tempeareture.
- ✓ 6. Mark your Samples with the Liquid Blocker Pen by drawing circles around them
- ✓ 7. Incubate with 10mM ethanolamine for 5 min at RT, if needed
- ✓ 8. Wash 3X 1 min with PBS
- ✓ 9. Incubate with 0.1% Triton-X100 in PBS for 1 min at RT
- ✓ 10. Wash 3X 1 min with PBS
- ✓ 11. Incubate the Samples with 1% NABH₄ in PBS , to deminish backgroundby dropping into the circles for 20 minutes (perhaps Liquid Blocker Pen circles need to be drawn again) **(only for tissue)**
- ✓ 12. Wash the samples 3 times with IF-buffer for 5 minutes
- ✓ 13. Block the autoactivity of the Tissue with Aurion 1% Aurion BSA, 1% Serum (Depending on source of secondary antibody) (Blocking Buffer) in IF-bufferfor 1 hour
- ✓ 14. Dilute multiple 1st Antibody (usualy 1:100) in Blocking Buffer and incubate overnight at 4°C

2nd Day

- ✓ 15. Wash the Slides 3 times with IF-buffer for 15 minutes
- ✓ 16. Incubate with the correspondinf 2nd Antibodies diluted in Blocking Buffer for 1 hour at RT
- ✓ 17. Wash the Slides 3 times with IF-buffer for 15 minutes in dark.

- ✓ 18. Incubate 5 min with 2X DAPI in PBS in dark.
- ✓ 19. Wash the Slides 1 time with IF-buffer for 5 minutes in dark.
- ✓ 20. Mount the coverslip using Fluorescence mounting media.