

# SMPD1 activity measurement

---

## Introduction

Get started by giving your protocol a name and editing this introduction.

## Materials

### > Substrate buffer

- > 0,1M Sodium acetate buffer pH 5,2 (120µl glacial acetic acid + 79ml 0,1N sodium acetate to 100ml)
- > 0,2%(w/v) synthetic sodium taurocholate
- > 0,02% sodium azide

### > HMU-PC substrate

- > 1.32 mM 6-hexadecanoylamino-4-methylumbelliferyl-phosphorylcholine (HMU-PC) in substrate buffer
- > (Mol. Wt: 594.72) [10mg in 12.72ml\*heat for short duration at 60 deg-C until clear. Aliquots (1ml each) and store at -80 deg-C]

### > Stop buffer

- > 0.2 M glycine–NaOH buffer (pH 10.7) (910µl 5M NaOH + 25ml of 0.2M glycine~volume to 100ml)
- > 0.2% (w/v) sodiumdodecyl sulphate
- > 0.2% (w/v) Triton X-100

## Procedure

1. mix 10µl of homogenate with (standardized quantity of protein 10 µg for fibroblasts) with 20µl of HMU-PC substrate in substrate buffer (1:1 ratio)
2. Incubate for 1 hour
3. Measure fluorescence kinetics for 30 min (**optional!**) at Ex. 360 Em. 460nm
4. Add 200µl stop buffer
5. Measure fluorescence (**Endpoint!**)