

IHC Multiplex

Introduction

IHC multiplex protocol for immune and cancer cells markers plus protein of interest. paraffin-embedded cuts on Tomo slides were used.

Materials

- › Xylene
- › Ethanol
- › ddH₂O
- › Hematoxylin (T865.2)
- › Eosin (HT110216-500ML)
- › Aqueous mounting medium
 - › 70% Glycerol
 - › 30% ddH₂O
- › PBS
- › PBS-T
 - › 0.1% Tween in PBS
- › 10X Antigen Retrieval Buffer
 - › 100mM citric acid
 - › 20mM EDTA
 - › in 800ml ddH₂O
 - › pH to 6.2 with NaOH
 - › 0.5% Tween20
 - › fill up to total 1000ml (can be stored at 4°C)
- › 3% H₂O₂ (1A8Y:1)
- › Aurion BSA 10% (900.022)
- › Antibody stripping buffer
 - › 1% SDS
 - › 25mM Glycine
 - › set pH to 2
 - › in ddH₂O
- › Picro-sirius red stain solution
 - › 0.1% direct red 80 (365548-5G)
 - › 0.1% fast green FCF (F7252-5G)
 - › Dissolved in saturated picric acid solution 1.3% in H₂O (P6744-1GA)
- › Secondary Antibodies ready to use

› DAKO EnVision+/HRP, Rabbit (K4003)

› DAKO EnVision+/HRP, Mouse (K4001)

› Primary Antibodies and the sequence of staining: please look at the end of protocol

› ImmPACT AMEC (SK-4285)

› Fast Red Substrate Kit (ab64254)

› VectaMount Permanent Mounting Medium (H-5000)

Procedure

Day 1

Day 1						
	A	B	C	D	E	F
1	Deparaffinize	1x	10min	rt.	Xylene	
2		1x	10min	rt.	Xylene	
3		1x	10min	rt.	Xylene	
4	Rehydrate	2x	5min	rt.	100% Ethanol	
5		1x	5min	rt.	95% Ethanol	
6		1x	5min	rt.	70% Ethanol	
7		1x	5min	rt.	50% Ethanol	
8	H&E	1x	5min	rt.	ddH ₂ O	
9		1x	1min	rt.	Hematoxylin	
10			2min	rt.	Running tap water	
11		1x	4min	rt.	Eosin	
12		2x	10sec	rt.	ddH ₂ O	
13	Mounting				Aqueous mounting medium	
14	Scanning					
15	Remove coverslip	1x	>20min	rt.	PBS	Do not pull the coverslip off!
16	washing	1x	2min	rt.	ddH ₂ O	
17	Antigen Retrieval	1x	35min	125°C	1X Antigen Retrieval Buffer (ARB)	in a pressure cooker
18	Cool down		30min		Remaining in the ARB	
19	Washing		5min		ddH ₂ O	
20	All of the following steps are carried out in a humidified chamber in the dark					
21	H ₂ O ₂ Blocking	1x	20min	rt.	3% H ₂ O ₂	in the dark
22	Washing	1x	5min	rt.	ddH ₂ O	
23		1x	2min	rt.	PBS	
24	Blocking	1x	60min	rt.	Blocking buffer	see primary ABs table
25	Primary Ab		o.n.	4°C	1:x	see primary ABs table

Day2

Day 2						
	A	B	C	D	E	F
1	Washing	3x	7min	rt.	PBS-T	
2	Secondary Ab	1x	60min	rt.	Ready-to-use	Based on host of the primary Ab
3	Washing	3x	7min	rt.	PBS-T	
4	Staining development	1x		rt.	AMEC staining solution	Under the microscope
5	Stopping HPO reaction	1x	>1min	rt.	ddH2O	
6	Hematoxylin	1x	1min	rt.	hematoxylin	
7	Washing		2min	rt.	Running tap water	
8	Mounting				Aqueous mounting medium	
9	Scanning					
10	Remove coverslip	1x	>20min	rt.	PBS	Do not pull the coverslip off!
11	washing	1x	2min	rt.	ddH2O	
12	AMEC destaining	1x	1min	rt.	70% Ethanol	
13		1x	3min	rt.	95% Ethanol	
14		1x	1min	rt.	70% Ethanol	
15		1x	2min	rt.	ddH2O	
16	Antibody stripping	1-2x	20min	55°C	Antibody stripping buffer pH 2	in shaking waterbath (90 rpm), also see primary ABs table
17	PH balancing	1x	7min	rt.	Antibody stripping buffer before pH adjusting (pH 6-7)	
18	Washing	1x	5min	rt.	ddH2O	
19		1x	2min	rt.	PBS	
20	Blocking	1x	60min	rt.	Blocking buffer	see primary ABs table
21	Primary Ab		o.n.	4°C	1:x	see primary ABs table
22	Go back to step 1 in day 2 (washing) and continue the cycle of multiplex					

Primary Antibodies and the sequence of staining									
	A	B	C	D	E	F	G	H	I
1	Antibody	Reactivity	Host	Application	Company	#Cat	Blocking in Aurion BSA	Dilution in Aurion BSA	Stripping
2	SMPD1	H	R	Ihc-multiplex	LSBio	LS-C334919	5% 1h R.T.	1:500 in 3% o.n.	2x 20min
3	CD8	H	M	Ihc-multiplex	Cell Marque	108M-94	1% 1h R.T.	1:100 in 1% o.n.	1x 20min
4	CD206	H	M	Ihc-multiplex	R&D	MAB25341	3% 1h R.T.	1:50 in 3% o.n.	2x 20min
5	PHH3	H	R	Ihc-multiplex	cell signaling	cs: 53348	1% 1h R.T.	1:100 in 1% o.n.	1x 20min
6	CD68	H	M	Ihc-multiplex	Dako	M0876	2% 1h R.T.	1:100 in 2% o.n.	1x 20min
7	CK19	H	M	Ihc-multiplex	Dako	M0888	2% 1h R.T.	1:85 in 0.5% o.n.	stained with fast red no stripping
8	CD3	H	M	Ihc-multiplex	Dako	M7254	1% 1h R.T.	1:70 in 1% o.n.	1x 20min
9	CD45	H	M	Ihc-multiplex	Dako	M0701	1% 1h R.T.	1:70 in 1% o.n.	no stripping, only destaining of Amec
10	Picro-sirius red staining								

Picro-sirius red staining						
	A	B	C	D	E	F
1	Picro-sirius red staining	1x	60min	rt.	Picro-sirius red solution	in the dark
2	Washing	2x	30sec	rt.	ddH2O	
3	Dehydration	1x	30sec	rt.	50% Ethanol	
4		1x	30sec	rt.	70% Ethanol	
5		1x	30sec	rt.	95% Ethanol	
6		2x	30sec	rt.	100% Ethanol	
7		3x	5min	rt.	Xylene	
8	Mounting				VectaMount Permanent Mounting Medium	
9	Scanning					

Notes

1. CK19 AB is hard to strip, so we stained it with the Fast red (permanent stain), which has a different color from ImmPACT AMEC (washable stain). the following ABs were in different cell compartment.
2. After the last AB cycle, we only destain AMEC and procende with Picro-sirius red staining.